

**KAYIT**

SAAT :

Numara

STERİLİZATÖR

1

2

3

4

5

PROGRAM

1

Bowie&amp;Dick 134 C 3,5 dk.

7

Inst.Containers 135 C 1

2

Instruments 134 C 10 dk.

8

Therm.Sensitive Inst/Kau

3

Textile 134 C 8 dk.

9

Fragile Packaging Texti

4

Rubber/Kauçuk 125 C 20 dk.

10

Flash 134 C 5 dk.

5

ATNC Inst . 134 C 18 dk.

11

Vacuum Leak Test 134

6

ATNC textile 134 C 18 dk.

STERİLİZATÖR GÖREVLİSİ :

**KONTROL**

GEÇERLİ

HATALI

Kor

NEM KONTROL

İNDİKATÖR KONTROL

GRAFİK KONTROL

STERİL DEPO GÖREVLİSİ :

ONAYLAYAN :



**Demo Hastanesi**  
MERKEZİ BUHAR STERİLİZASYON KAYIT VE KONTROL  
FORMU



Doküman Kodu: DS.FR.005

Yayın Tarihi: 15.09.2021

Revizyon Tarihi: -

Revizyon No: 0

Sayfa: 2 / 2

(P 11 için yalnız grafik kontrolü yapılır)